

## Animal Islet Isolation Protocol


### **Materials Inventory and Preparation Instructions**

#### Surgical Instrument and Digestion Apparatus Pack

<b>Item</b>	<b>Quantity</b>	<b>Check</b>
Nalgene PP bin	1	
Medium SS tray	1	
Small SS tray w/lid	1	
Hammer SS	1	
Metzenbaum scissors blunt	1	
Metzenbaum scissors sharp	1	
Gemini clamp	1	
Mayo scissors	1	
Classic Crile forceps	2	
Delicate Crile forceps straight	3	
Delicate Crile forceps curved	4	
Dressing forceps	2	
Children's dressing forceps	1	
Ricordi chamber (250 ml)	1	
SS screen 381 um	1	
SS screen 450 um	1	
Silicon nitride ½" balls	7	
Heating coil	1	
100 ml wide mouth graduate	1	
1 L Erlenmyer flasks	2	
Surgical gauze in a sterilization pouch	1	
Digestion circuit L/S 16	1	
Digestion circuit L/S 17	1	
Thermometer clamps	4	
Three prong fixed clamps	3	
Three prong adjustable clamp	1	

## Preparation of Surgical Instruments and Digestion Apparatus

- 1) Wipe down the counter with IPA.
- 2) Place two layers of surgical wrapping on the counter.
- 3) Place the autoclave tray (or PP bin) on the wrapping on the counter.
- 4) Place the stainless steel pans inside the autoclave tray (or PP bin).
- 5) Place the surgical instruments inside the instrument tray, cover, and place in the autoclave tray.
- 6) Place the Ricordi digestion chamber with screen and 5 marbles, stand, and heating coil in the autoclave tray.
- 7) Place 20 pieces of 4"x4" gauze in an autoclave pouch, seal, and place in the autoclave tray.
- 8) Aluminum foil wrap the top of a 100 ml wide-mouth graduated cylinders and place in the autoclave tray.
- 9) Aluminum foil wrap the top of one 1 L Erlenmeyer flask and one 400 ml Pyrex beaker and place in the autoclave tray.
- 10) Cut the L/S 16 tubing into the following sizes (2 x 4.5 in, 1 x 11.5 in, 1 x 10 in, 1 x 20 in, 1 x 50 in, 1 x 65.0 in). Following the diagrams shown below, connect tubing to the Y and T connectors, place pinch clamps, place cable ties on all connections, and place tubing set into a sterilization pouch; seal and label the pouch and place in the autoclave tray.
- 11) Cut the L/S 17 tubing into the following sizes (1 x 5 in, 1 x 11.0 in, 1 x 26 in, 1 x 42.0 in). Following the diagrams shown below, connect tubing to the 3-way stopcock and T-connector, place cable ties on all connections and place into a sterilization pouch; seal and label the pouch and place in the autoclave tray.
- 12) Place three fixed 3-prong ring clamps, 1 adjustable 3-prong clamp, and 4 thermometer clamps in the autoclave tray.
- 13) Place a cleaned media collection bell and tubing adaptor in an autoclave pouch, seal, and place in the autoclave tray.
- 14) Wrap the autoclave tray as for a sterile surgical instrument pack.
- 15) Tape the outside of the package with autoclave indicator tape and label with contents, date, and name.
- 16) Sterilize the surgical instrument pack by autoclave using a PreVac cycle at 121<sup>0</sup>C for 20-30 minutes.
- 17) Store the complete/sealed Bin in room G24A.



## **Disposables Pack**

<b>Item</b>	<b>Quantity</b>	<b>Check</b>
Disposable scalpel	1	
Monotherm temp probe	4	
Monotherm temp cables	2	
2.0 Suture	1	
3 Way stopcock	1	
5 cc syringe	2	
10 cc syringe	2	
60 cc syringe	4	
18 ga Angiocath	2	
20 ga Angiocath	2	
24 ga Angiocath	2	
18 ga Needle	5	
20 ga Needle	5	
Petri dishes 10 x 35	1 sleeve	
Petri dishes 10 x 35 w/grid	1 sleeve	
Petri dishes 100 mm	1 sleeve	
Tissue culture plate, 12 well	1 pack	

## **Islet Purification Pack**

<b>Item</b>	<b>Quantity</b>	<b>Check</b>
Classic Crile forceps straight (on Cobe)	5	
Blood processing bags	2	
COBE couplers	2	
60 cc syringe	2	
Fenwal Transfer bag 600 ml	2	

## Reagents

<b>Item</b>	<b>Quantity</b>	<b>Check</b>
2 % FBS in HBSS (Flushing solution)	1 L	
2.5% FBS in UW (Islet suspension/wash soln.)	1 L	
2.5% FBS in RPMI (Dilution solution)	5 L	
CMRL 1066 w/10% FBS (Islet culture)	1 L	
1x HBSS	1 L	
Decontamination soln 1		
Decontamination soln 2		
Decontamination soln 3		
Biocoll BL (130 ml)	1	
Biocoll CL 1 (130 ml)	1	
Biocoll CL 2 (140 ml)	1	
Krebs buffer (static incubation)	1 L	
FDA stock	1	
PI stock	1	
DPBS	1 L	
Dithizone soln	50 ml	
1x PBS sheath buffer	25 L	
Euro-Ficoll 1.037	75 ml	
Euro-Ficoll 1.096	75 ml	
Euro-Ficoll 1.108	75 ml	

## Complete directions for media preparation can be found in UW03 ISL-SOP-002.

## Islet Isolation Procedure

1. Adult Beagle or Rhesus – surgically euthanized with pancreas perfused *in situ* with cold UW solution.
2. Pancreas is stored in 250 ml of ice-cold UW solution and transported directly to the islet isolation lab.
3. Islets are isolated using a modified version of the Miami Diabetes Research Institute protocol for canine islet isolation.
  - a. Pancreas cannulated and then distended and digested with Liberase CI in 1x HBSS (Canine [0.25 g], Rhesus [0.125 g] /500 ml total vol.) in a Ricrodi digestion chamber and recirculation circuit at 37°C for up to 30 minutes.
  - b. Monitor digest progression by taking samples, staining with dithizone solution and observing with the inverted microscope.
  - c. Dilution of pancreas digest is accomplished with 6-7 L of 10% RPMI and collection is directly into 500 ml conical centrifuge tubes containing 30 ml of 2.5% FBS in RPMI.
  - d. Digest is pelleted by centrifugation at 1000 rpm for 2 minutes at 4°C.
  - e. After the spin, the supernatants are poured into a waste beaker or aspirated into the waste flask. The pellets are gently resuspended and combined into a separate 500 ml conical tube, which is stored on ice until all of the collected digest has been centrifuged.
  - f. The final collected digest is washed in 2.5% FBS in UW solution and pelleted by centrifugation at 1000 rpm for 2 minutes at 4°C.
  - g. The supernatant is discarded and the pellet is resuspended in 2.5% FBS in UW solution according to the packed tissue volume.
    - i. 15 – 20 ml of tissue – bring to 100 ml with 2.5% FBS in UW and run gradient on the COBE.
    - ii. 30 – 50 ml of tissue – bring to 200 ml with 2.5% in UW and run two gradients on the COBE.
  - h. Refer to UW03 ISL-CP for gradient loading and COBE centrifugation instructions.
    - i. After 2-3 minutes of centrifugation at 2000 RPM collect 17, 30 ml fractions into 50 ml conical centrifuge tubes (DO NOT STOP THE COBE).
    - j. Take 100 µl samples of each fraction, transfer to a 10 x 35 mm counting dish containing a small quantity of PBS, stain with dithizone, and observe with the inverted microscope for islets.
    - k. Determine islet purity, and yield in IEQ manually. Data are recorded on the Data Summary Sheet.

## **Modified Islet Purification Method – Discontinuous Gradient Centrifugation**

1. Pancreas digest is washed in Eurocollins solution.
2. Digest is resuspended in 1.108 g/ml Euro-Ficoll gradient at a ratio of 1:11 (tissue: Ficoll) and 12 mL is aliquoted per 50 ml conical centrifuge tubes.
3. Ten mL of 1.096 g/ml and ten ml of 1.037 g/ml Euro-Ficoll are then added to each discontinuous gradient.
4. The gradients are centrifuged at 900 xg for 12 – 15 minutes and the islets collected at the interfaces.
5. Purified islets are washed in CMRL culture media, pooled, and counted.

10 mL 1.037 g/mL EF
10 mL 1.096 g/mL EF
12 mL Islets in 1.108 g/ml EF